Manufactured By: HYPHEN BioMed

# Human Thrombin Thrombine humaine

Vial of 100 NIH/flacon de 100 NIH  
Set of 10x100 NIH / Coffret de 10x100 NIH

Ref. AEZ006A/AEZ006L

STORE AT  
2-8 °C

For research use only

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CONSERVER A  
2-8 °C

**Origin:**

Prepared from a prothrombin complex concentrate extracted from human citrated plasma and fully activated in a purified system. Thrombin is then purified using an ion exchange chromatography and stabilised. It is a highly purified preparation, mainly in the  $\alpha$  form, with a specific activity of more than 2,000 NIH/mg.

**Presentation:**

Thrombin activity is reported in harmonized NIH or IU or USP units by reference to the WHO/NIBSC International Standard for Human Thrombin (01/580) (refer to "A reunification of the US ("NIH") and International Unit into a single standard for Thrombin", Longstaff et al, Thromb Haemost 2005, 93:261-6)). NIH is a clotting unit for thrombin activity. Stabilised and lyophilised human Thrombin. Vial containing approximately 100 NIH of thrombin. AEZ006A: 1x100 NIH ; AEZ006L: 10x100 NIH

**Reconstitution:**

Each vial must be restored with 1 ml distilled water, in order to get a concentration of about 100 NIH/ml. Dilute in the adequate buffer, if necessary.

**Excipients:** Content per vial: Poly-Ethylene-Glycol 6000 (PEG 6000) (20mg), Sodium Chloride (6mg), Glycine (20mg), BSA (10mg), stabilizers.

**Purity:** Human Thrombin has one major band of 35,000 daltons on SDS-PAGE.

**Activity:**

- Clotting time (at 10NIH/ml) on purified human Fibrinogen (5 mg/ml):  $5'' \pm 1''$
  - Clotting time (at 10NIH/ml) on citrated human plasma pool:  $5'' \pm 1''$
  - Chromogenic activity on CS-01(38) substrate:  $\geq 0.75/\text{min.NIH.ml}$
  - Clotting activity determined in NIH (or IU).
  - Chromogenic activity on specific substrates (expressed in nkats):
- When tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the following specific activities are observed on the various IIa substrates (expressed in nkats/ $\mu\text{g}$ ):

CS-01(38)	CS-01(81)
2.9 $\pm$ 1.0	3.6 $\pm$ 1.0

The exact activity in nkats/ $\mu\text{g}$  with CS-01(38) is reported on the analysis certificate for each lot.

**Viral safety:** The human plasma used for Thrombin purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. Bovine Serum Albumin (BSA) was prepared from bovine plasma, which was tested for the absence of infectious agents, and collected from animals free from BSE. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

**Use:** For *in vitro* use only. All research studies and protocols where a source of highly purified human Thrombin is necessary.

**Stability:**

- Restored :**
- 7 days at room temperature (18-25 °C).
  - 21 days at 2-8 °C.
  - 6 months at -30 °C or below.
- Lyophilized :** Until the expiration date printed on the vial.

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