Manufactured By: HYPHEN BioMed

Purified human Factor Xa
Facteur Xa humain purifié
Vial of 0.1 mg / flacon de 0.1 mg
Set of 6x0.1 mg / coffret de 6x0.1 mg

Ref. AEZ007B / AEZ007L

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified factor X from citrated human plasma, activated with agarose bound RVV. MW of about 44,000.

Presentation:

Vial containing approximately 0.1 mg of highly purified human factor Xa in presence of additives and preservatives. (1 Unit FXa is the amount which can be generated from 1 ml citrated human plasma, ie about 10 µg).

Reconstitution:

Each vial must be restored with 1 ml distilled water, in order to get a concentration of about 100µg/ml. It can be diluted to the desired concentration, when required, in the adequate buffer.

Excipients:

Content per vial: Glycine (20mg), Hepes (6mg), Prionex (0.5% after reconstitution), PEG6000 (5 mg), Sodium chloride (9mg).

Purity:

The factor X zymogen used for preparing factor Xa has one major band of 59,000 daltons on SDS-PAGE.

Activity:

a) Tested for its amidolytic activity on a FXa specific substrate, in a purified system without or with RVV. All FX is converted to FXa.
b) Chromogenic activity on specific substrates (expressed in nkats):
When tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the following specific activities are observed on the various FXa substrates (expressed in nkats/µg):

Table with 3 columns: CS-11(22), CS-11(32), CS-11(65) and 3 rows of activity values.

The exact activity in nkats/µg with CS-11(32) is reported on the analysis certificate for each lot.

Viral safety:

The human plasma used for Factor Xa preparation was tested with registered methods and found negative for HIV (1 and 2) antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only. All research studies and protocols where a source of highly purified human Factor Xa is necessary. Studies on heparin and Low Molecular weight heparin.

Stability:

- Restored: - 7 days at room temperature (18-25 °C)
- 1 month at 2-8 °C
- 6 months at -20 °C or below.

Lyophilized: Until the expiration date printed on the vial.

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