Manufactured By: HYPHEN BioMed

Purified Human Factor XIIa
Facteur XIIa humain purifié
Vial of 100 ng (AEZ012A) or Set of 6x100 ng (AEZ012K)
Flacon de 100 ng (AEZ012A) ou Coffret de 6x100 ng (AEZ012K)

Ref. AEZ012A / AEZ012K

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified Factor (α) XIIa (also known as Prekallikrein Activator) from citrated human plasma, activated from human FXII using Dextran sulphate, removed after activation.

Presentation:

Vial containing approximately 100 ng (i.e. about 80 IU*) of highly purified human factor XIIa in presence of additives and preservatives.

(*International Unit according to the 2nd NIBSC International Standard for Prekallikrein Activator code O2/168).

Reconstitution:

Each vial must be restored with 1 ml distilled water. It can be diluted to the desired concentration, when required, in the adequate buffer.

Excipients:

- Tris, Bovine Serum Albumin (BSA).
- 9 mg sodium chloride/vial.

Purity:

Purity of the Factor XII zymogen (one major band expected at about 80,000 daltons) used for preparing Factor XIIa is controlled on SDS-PAGE, as well as complete activation.

Activity:

Tested for its amidolytic activity in a FXIIa chromogenic assay.
100ng correspond to about 80IU.

Viral safety:

The human plasma used for Factor XIIa preparation was tested with registered methods and found negative for HIV (1 and 2) antibodies, HBs Ag and HVC antibodies. BSA was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.
All research studies and protocols where a source of highly purified human Factor XIIa is necessary.

Stability:

- Restored :**
- 24 hours at room temperature (18-25 °C)
 - 72 hours at 2-8 °C
 - 2 months at -20 °C or below

Lyophilized: Until the expiration date printed on the vial.

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