Manufactured By: HYPHEN BioMed

Purified Human Fibrinogen, Plasminogen Free
Fibrinogène humain purifié sans Plasminogène
Vial of 5 mg (APP002A) /flacon de 5 mg (APP002A)
Set of 6x5 mg (APP002K) /Coffret de 6x5 mg (APP002K)

Ref. APP002A / APP002K

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified fibrinogen extracted from citrated human plasma. Plasminogen is specifically removed by affinity chromatography.

Presentation:

Vial containing approximately 5 mg of highly purified human fibrinogen, plasminogen free.

Reconstitution:

Each vial must be restored with 1 ml distilled water in order to obtain a 5 mg/ml fibrinogen solution. This solution can be diluted to the desired concentration when required.

Excipients:

- Glycine, trisodium citrate and sodium veronal buffer
- 6 mg sodium chloride/vial

Purity:

One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clottability $\geq 98\%$.

Activity:

Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/ml (100 μ l), is clotted by thrombin, at 5 NIH/ml (100 μ l), the clotting time must be below 15 sec.

Absence of plasminogen is tested with a clot lysis method in presence of tPA

Viral safety:

The plasma used for fibrinogen purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.

All research studies and protocols where a source of highly purified fibrinogen is necessary.

Stability:

- Restored :**
- 48 hours at room temperature
 - 96 hours at 2-8 °C

Lyophilized : Until the expiration date printed on the vial.

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Vial of 5 mg (APP002A) /flacon de 5 mg (APP002A)
Set of 6x5 mg (APP002K) /Coffret de 6x5 mg (APP002K)

Ref. APP002A / APP002K

STORE AT
2-8 °C

Lot:

Expiration:

CONSERVER A
2-8 °C

ANALYSIS CERTIFICATE

ANALYTICAL DATA	SPECIFICATIONS
1. Protein Content (Lowry method; A280 nm) mg/vial	> 4.5mg/vial
2. SDS-PAGE (5 % acrylamide) 1 single band of about daltons	1 major band of about 340,000 daltons
3. Clotting time (5 NIH/ml Thrombin) sec	< 15 sec
4. Coagulability %	≥ 98%
5. Clot lysis in presence of tPA No lysis	No lysis

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DATE :

CONCLUSIONS : Passed/Conforme Refused/Non-conforme

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