Manufactured By: HYPHEN BioMed

Anti-Thrombin (AT) Anti-Thrombine (AT)

Vial of 150 µg (APP004A) / Flacon de 150 µg (APP004A)
Set of 6x150 µg (APP004K) / Coffret de 6x150 µg (APP004K)

Ref. APP004A / APP004K

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified from human plasma using a combination of affinity and ion-exchange chromatographies.

Presentation:

Vial containing approximately 150 µg of AT corresponding to about 1 Plasma Equivalent Unit or PEU (amount of AT in 1 ml human plasma), which is about 1 International Unit (IU) (when determined and standardized against the WHO/NIBSC International standard for Antithrombin, plasma, human).

Reconstitution:

Each vial must be restored with 1 ml distilled water, in order to get a concentration of about 150 µg/ml. Dilute in the adequate buffer, if necessary.

Excipients:

Content per vial: Glycine (20mg), Hepes (6mg), Sodium Chloride (9mg).

Purity:

One major band at about 58,000 daltons on SDS-PAGE

Activity:

In a clotting assay for the thrombin inhibition capacity in presence of an excess of unfractionated heparin : 150 µg AT neutralize at least 200 NIH units of thrombin. Tested for anti-Xa activity in presence of heparin, with a chromogenic assay (> 6IU/mg).

Factor Xa neutralisation capacity is > 500 µg of Xa/mg AT.

The absence of heparin is verified after purification, using a toluidine blue assay (refer to: F.C MacIntosh, "A colorimetric method for the standardization of heparin preparations", Biochem, 1941, 35:776-782), and an anti-Xa specific chromogenic assay.

The exact activity is reported on the analysis certificate for each lot.

Viral safety:

The plasma used for AT purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For *in vitro* use only. All research studies and protocols where a source of highly purified human AT is necessary.

Stability:

Restored :

- 24 hours at room temperature (18-25 °C)
- 72 hours at 2-8 °C
- 2 months frozen at -30 °C or below

Lyophilized : Until the expiration date printed on the vial.

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