Manufactured By: HYPHEN BioMed

Fibrin Fragment DD

Fragment DD de la Fibrine

Vial of 200 µg/flacon de 200 µg

Ref. APP017A

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified fibrin Fragment DD, obtained by degradation of a Fibrin clot with plasmin, then purification with ion exchange and gel filtration chromatographies.

Presentation:

Vial containing approximately 200 µg of highly purified human fibrin, fragment DD (DDimer).

Reconstitution:

Each vial must be restored with 1 ml distilled water. It can be diluted to the desired concentration when required. Use a buffer containing a carrier protein (such as BSA) for dilution, in order to avoid material loss.

Excipients:

- Glycine, Hepes.
- 9 mg sodium chloride/vial

Purity:

One major band of about 160,000 daltons on SDS-PAGE.

Viral safety:

Fibrinogen used for preparing fragment DD was extracted from human plasma tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.
All research studies and protocols where a source of highly purified fibrinogen is necessary.

Stability:

- Restored :**
- 72 hours at room temperature
 - 7 days at 2-8 °C
 - 6 months at -20 °C or below

Lyophilized : Until the expiration date printed on the vial.

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