Manufactured By: HYPHEN BioMed

Fibrinogen Fragment E3t

Fragment E3t du Fibrinogène

Vial of 200 µg/flacon de 200 µg

Ref. APP019A

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Highly purified fibrinogen Fragment E3t (terminal degradation stage), obtained by extensive degradation of Fibrinogen with plasmin, then purification with ion exchange and gel filtration chromatographies.

Presentation:

Vial containing approximately 200 µg of highly purified human fibrinogen, fragment E3t (terminal degradation stage).

Reconstitution:

Each vial must be restored with 1 ml distilled water. It can be diluted to the desired concentration when required. Use a buffer containing a carrier protein (such as BSA) for dilution, in order to avoid material loss.

Excipients:

- Glycine, Hepes.
- 9 mg sodium chloride/vial

Purity:

One major band of about 45,000 daltons on SDS-PAGE.

Viral safety:

Fibrinogen used for preparing fragment E3t was extracted from human plasma, tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.
All research studies and protocols where a source of highly purified fibrinogen fragment E3t is necessary.

Fibrinogen degradation produces first fragment X, then fragments Y and D and finally fragments D and E. With progression of degradation first fragment E1 is produced, then fragment E2 and finally fragments E3 and E3t.

Stability:

- Restored :
- 72 hours at room temperature
 - 7 days at 2-8 °C
 - 6 months frozen at -20 °C or below

Lyophilized : Until the expiration date printed on the vial.

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