Manufactured By: HYPHEN BioMed

Recombinant rat PAI-1 PAI-1 de rat recombinant

Vial of 10 µg/flacon de 10 µg

Ref. ARE0050

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Recombinant rat PAI-1 expressed in E. Coli.

Presentation:

Vial containing 10 µg of recombinant rat PAI-1 (expressed in E. Coli), lyophilized.

Reconstitution:

Restore the vial with 1 ml of distilled water.

Nota : If PAI-1 must be used at a concentration below 0.01 mg/ml, the dilution buffer must contain Bovine Serum Albumin (i.e. 1%) or any other carrier protein or polymer in order to avoid loss of the protein by adsorption on the container walls.

Excipients:

– Sodium Phosphate, Arginine, Bovine Serum Albumin, Sodium Chloride, pH 7.50

Activity:

Rat-PAI-1 has an activity of about 50% (tested for its binding to tPA). PAI-1 activity is labile. Rat-PAI-1 must be used rapidly or stored frozen in small aliquots and thawed just before use.

Characteristics:

PAI-1 has a molecular weight of about 50,000 daltons.

The Rat-PAI-1 preparation has a purity grade > 95%.

Use:

For in vitro use only.
All research studies and protocols where a source of recombinant rat PAI-1 is necessary.

Viral safety:

The bovine plasma used for BSA preparation was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

Stability:

Restored :

- 4 hours at room temperature
- 24 hours at 2-8 °C
- 6 months at -20 °C or below

Lyophilized : - Until the expiration date printed on the vial.

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