# PLASMA PAI-1 CONTROL

## # ASC011K

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Last revision / Dernière version: 12/01/2006

### Intended use:

Lyophilised plasmas, at defined PAI-1: Ag and PAI-1 Activity concentrations, as measured with the ZYMUTEST-PAI-1: Ag kit (ref ARK012A), or the Zymutest PAI-1 Activity kit (ref ARK019A) for the quality control of PAI-1 assays.

### Reagents:

Each kit contains:

- 3 vials of 1 ml of normal plasma, lyophilised, at a high PAI-1 concentration (Control I)
- 3 vials of 1 ml of normal plasma, lyophilised, at a low PAI-1 concentration (Control II).

The PAI-1: Ag and PAI-1 Activity concentrations and the acceptancy ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

### Reagent Preparation and Stability:

In their original package, and before any use, when stored at 2-8 °C, the reagents are stable until the expiration date printed on the kit.

### Excipients:

Glycine, stabilizers.

### Preparation:

- **PAI-1 Control I (high):** restore the vial with 1 ml distilled water; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **PAI-1 Control II (low):** restore the vial with 1 ml distilled water; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These plasmas are stable:

### When used for Antigenic assays:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

### When used for Activity assays:

- 1 hour at room temperature (18-25 °C)
- 4 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

**Nota :** Plasmas used for PAI-1 controls were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious. (win HIV/HBS).

### Application

The PAI-1: Ag, the PAI-1 Activity concentrations, indicated on the flyer, for Control I and Control II, have been measured with the ZYMUTEST PAI-1: Ag (ref ARK012A) or the Zymutest PAI-1 Activity (ref ARK019A) kits.

These assays are standardised by reference to the NIBSC international standard.

The target values and the acceptancy ranges, indicated on the flyer, are valid for ZYMUTEST PAI-1: Ag and Zymutest PAI-1 Activity. The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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