**PLASMA TAFI CONTROL**  
# ASC014K

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Last revision / Dernière version : 12/01/2006

**Intended use**

Lyophilised plasmas, at a defined TAFI: Ag concentration, as measured with the ZYMUTEST-TAFI Ag kit (ref ARK008A), for the quality control of TAFI antigenic assays.

**Reagents:**

Each kit contains:

- 3 vials of 0.5 ml of normal plasma, lyophilised, at a high TAFI concentration (Control I)
- 3 vials of 0.5 ml of normal plasma, lyophilised, at a low TAFI concentration (Control II).

The TAFI: Ag concentrations and the acceptancy ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

**Reagent Preparation and Stability:**

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

**Excipients:**

Glycine, stabilizers.

**Preparation:**

- **TAFI Control I (high):** restore the vial with 0.5 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **TAFI Control II (low):** restore the vial with 0.5 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These plasmas are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

**Nota :** Plasma used for TAFI controls were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious. (win HIV/HBS).

**Application**

The TAFI: Ag concentration, indicated on the flyer, for Control I and Control II, has been measured with the ZYMUTEST TAFI: Ag (ARK008A), Elisa kit.

The target values and the acceptancy ranges, indicated on the flyer, are valid for ZYMUTEST TAFI: Ag only. The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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