PLASMA UTA CONTROL

ASC017K

Manufactured By: HYPHEN BioMed

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Last revision / Dernière version: 22/10/2007

Intended use

Lyophilised plasmas, at defined Annexin V, tPA: Ag and uPA: Ag concentrations, as measured with the ZYMUTEST-Annexin V kit (ref ARK004A), Zymutest tPA: Ag kit (ref ARK011A), Zymutest uPA: Ag kit (ref ARK013A) for the quality control assays.

Reagents:

Each kit contains:

- 3 vials of 1 ml of normal plasma, lyophilised, at a high Annexin V, tPA and uPA concentration (Control I)
- 3 vials of 1 ml of normal plasma, lyophilised, at a low Annexin V, tPA and uPA concentration (Control II).

The concentrations and the acceptance ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

Reagent Preparation and Stability:

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

Excipients:

Glycine, stabilizers.

Preparation:

- **UTA Control I (high):** restore the vial with 1 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **UTA Control II (low):** restore the vial with 1 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These plasmas are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota : Plasmas used for UTA controls were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Application

The Annexin V, tPA: Ag and uPA:Ag concentrations, indicated on the flyer, for Control I and Control II, have been measured with the corresponding ZYMUTEST (ARK004A, ARK011A, ARK013A), Elisa kit.

This assay is standardised by reference to the NIBSC international standard.

The target values and the acceptance ranges, indicated on the flyer, are valid for ZYMUTEST Annexin V, tPA: Ag ,uPA:Ag only. The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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