HIGH POSITIVE DDIMER CONTROL PLASMA
ASC101K

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Last revision / Dernière version: 30/07/2009

Intended use

Lyophilised plasma, at a defined DDimer concentration, as measured with commercially available methods and especially the VIDAS-DDimer assay.

Reagents:

Each kit contains:

- **C1:** Vials of 1 ml of High Positive DDimer Control Human Plasma, lyophilised, at a positive DDimer concentration.

The DDimer concentration and the acceptance range, for this positive control, can present variations from lot to lot, but are exactly indicated on the flyer provided with each set of positive control.

The DDimer concentration measured with the Vidas DDimer kits is of about **2 µg/ml**.

Reagent Stability:

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

Excipients:

Glycine, stabilizers.

Preparation:

- **High Positive DDimer Control Human Plasma (C1):** restore the vial with 1 ml of distilled water; let for 15 min at room temperature; mix gently until complete dissolution of the content.

This positive DDimer control plasma is stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in its original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota : Plasmas used for preparing the Positive DDimer control plasma were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Application

The DDimer concentration, indicated on the flyer, for High Positive DDimer Control plasma, has been measured with the Vidas DDimer, Elisa (Biomerieux).

The target value and the acceptance range, indicated on the flyer, are valid for the assay used. The results obtained can vary, when this control is used with other kits, according to the assay reactivity and its standardisation. This control can nevertheless be used with different DDimer assays, for the quality control, and for testing run to run variability, but the acceptance value and range must be determined by each laboratory.

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