Imidazole buffer

REF	AR021A-RUO,	Bf	1x25 mL
REF	AR021K-RUO,	Bf	3x25 mL
REF	AR021L-RUO,	Bf	12x25 mL
REF	AR021B-RUO,	Bf	1x15 mL
REF	AR021M-RUO,	Bf	3x15 mL
REF	AR021N-RUO,	Bf	12x15 mL

Imidazole buffer
FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.

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English, last revision: 10-2020

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INTENDED USE:

Dilution buffer for coagulation assays.
Ready to use.
This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

Bf Solution containing Imidazole (0.05M), sodium chloride (0.1M) and small amounts of Sodium azide (0.2 g/L) as preservative, pH 7.30 ± 0.20. Ready to use.

REF	AR021A-RUO → 1 vial of 25 mL
REF	AR021K-RUO → 3 vials of 25 mL
REF	AR021L-RUO → 12 vials of 25 mL
REF	AR021B-RUO → 1 vial of 15 mL
REF	AR021M-RUO → 3 vials of 15 mL
REF	AR021N-RUO → 12 vials of 15 mL

WARNING AND PRECAUTIONS:

- In contact with lead or copper pipes, sodium azide can generate explosive compounds.
- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- This device of *in vitro* use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Bf Reagent is ready to use, homogenize and load it on the analyzer following application guide instruction.
For another use, allow to stabilize for 30 minutes at room temperature (18-25°C), homogenize before use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

Bf Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of:

- **8 weeks** at 2-8°C.
- **8 weeks** at room temperature (18-25°C).
- **Stability on board of the analyzer: see the specific application.**

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the test to be performed.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

H360D: May damage the unborn child.

Changes compared to the previous version.

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