Purified bovine Factor Xa

REF BE101D Ez Vial of 30 µg / Flacon de 30 µg

REF BE101L Ez 10 vials 30 µg / 10 flacons de 30 µg

Purified bovine Factor Xa
FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES

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English, revision: 09-2022

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INTENDED USE:

All research studies and protocols where a source of highly purified bovine Factor Xa (FXa) is necessary. Studies on Heparin and Low Molecular Weight Heparin.

This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

Ez Highly purified Factor X (MW is about 55KDa) from bovine plasma, activated with agarose bound RVV. FXa activity is expressed in unit (U), related to **WHO reference reagent for the FXa Human (h) in force** (NIBSC), lyophilized. Contains approximately 30 µg of highly purified bovine FXa, BSA and stabilizers.

The product is classified as non-hazardous and is not subject to labeling according to EC Regulation No. 1272/2008 [CLP].

WARNING AND PRECAUTIONS:

- This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.
- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- This device of *in vitro* use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Gently remove the freeze-drying stopper, to avoid any product loss when opening the vial.

Ez Reconstitute the contents of each vial with exactly **10 mL of distilled water** in order to get a FXa concentration of about 3 µg/mL. It can be diluted to the desired concentration, when required, in an adequate buffer such as Tris NaCl BSA 1% pH 7.40. Shake vigorously until complete dissolution while avoiding formation of foam and allow to stabilize for 15 minutes at room temperature (18-25°C), homogenize before use.

Activity:

Raw material is tested for its amidolytic activity on a FXa specific substrate, in a purified system without or with RVV. All FX is converted to FXa. (A405(RVV-) / A405 (RVV+) >0.90). Chromogenic activity of final product on CS11-(6S) substrate (expressed in U and nkats for information): tested in the optimized buffer conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl).

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

Ez Reagent stability after reconstitution, free from any contamination or evaporation, and stored closed, is of:

- 3 months** at 2-8°C
- 7 days** at room temperature (18-25°C)
- 6 months** frozen at -20°C or less* for all reference

*Thaw only once, as rapidly as possible at 37°C and use immediately.

LIMITATIONS:

- Any reagent presenting no limpid appearance or showing signs of contamination must be rejected. **The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

Changes compared to the previous version.

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