# Human Thrombin

## Thrombine humaine

Vial of 10 NIH/flacon de 10 NIH  
Set of 6x10 NIH / Coffret de 6x10 NIH

Ref. EZ0060 /EZ006K

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English, last revision 12-2016

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**ORIGIN:**

Prepared from a prothrombin complex concentrate extracted from human citrated plasma and fully activated in a purified system. Thrombin is then purified using an ion exchange chromatography and stabilised. It is a highly purified preparation, mainly in the  $\alpha$  form, with a specific activity of more than  $\geq 1,500$  NIH/mg.

**PRESENTATION:**

Thrombin activity is reported in harmonized NIH or IU or USP units by reference to the WHO/NIBSC International Standard for Human Thrombin (01/580) (refer to "A reunification of the US ("NIH") and International Unit into a single standard for Thrombin", Longstaff *et al*, Thromb Haemost 2005, 93:261-6). NIH is a clotting unit for thrombin activity. Stabilised and lyophilised human Thrombin.

Vial containing approximately 10 NIH of thrombin.

EZ0060: 1x10 NIH ; EZ006K: 6x10 NIH.

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**RECONSTITUTION:**

Each vial must be restored with 1 mL distilled water, in order to get a concentration of about 10 NIH/mL. Dilute in the adequate buffer, if necessary.

**EXCIPIENTS:**

Content per vial: Poly-Ethylene-Glycol 6000 (PEG 6000) (20mg), Sodium Chloride (6mg), Glycine (20mg), BSA (10mg), Calcium, stabilizers.

**PURITY:**

Human Thrombin has one major band of 35 000 daltons on SDS-PAGE.

**ACTIVITY:**

Clotting time on purified human Fibrinogen (4 mg/mL):  $8'' \pm 3.5''$ .

Clotting time on citrated human plasma pool:  $5'' \pm 1''$ .

Chromogenic activity on CS-01(38) substrate:  $\geq 0.75/\text{min.NIH.mL}$ .

Clotting activity determined in NIH (or IU).

Chromogenic activity on specific IIa substrates (expressed in nkats/ $\mu\text{g}$ ) tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the exact activity in nkats/ $\mu\text{g}$  with CS-01(38) is reported on the analysis certificate for each lot.

**VIRAL SAFETY:**

The human plasma used for Thrombin purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. Bovine Serum Albumin (BSA) was prepared from bovine plasma, which was tested for the absence of infectious agents, and collected from animals free from BSE. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

**USE:**

For *in vitro* use only. All research studies and protocols where a source of highly purified human Thrombin is necessary.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

Restored : - 7 days at room temperature (18-25°C).

- 21 days at 2-8°C.

- 6 months at -30°C or below.

Lyophilized : Until the expiration date printed on the kit.

The stability of the reconstituted reagent should be checked under laboratory work conditions.

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