# Purified human Factor XIIa

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Vial of 100 ng (EZ012A)  
Flacon de 100 ng (EZ012A)

Ref. EZ012A

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English, last revision 08-2017

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**ORIGIN:**

Highly purified Factor ( $\alpha$ ) XIIa (also known as Prekallikrein Activator) from citrated human plasma, activated from human FXII using Dextran sulphate, removed after activation.

**PRESENTATION:**

Vial containing approximately 100 ng (i.e. about 80 IU\*) of highly purified human Factor XIIa in presence of additives and preservatives.

(\*International Unit according to the NIBSC International Standard in force).

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**RECONSTITUTION:**

Each vial must be restored with 1 mL distilled water. It can be diluted to the desired concentration, when required, in the adequate buffer.

**EXCIPIENTS:**

Tris, Bovine Serum Albumin (BSA).  
9 mg sodium chloride/vial.

**PURITY:**

Purity of the Factor XII zymogen (one major band expected at about 80,000 daltons) used for preparing Factor XIIa is controlled on SDS-PAGE, as well as complete activation.

**ACTIVITY:**

Tested for its amidolytic activity in a FXIIa chromogenic assay.  
100ng correspond to about 80IU.

**VIRAL SAFETY:**

The human plasma used for Factor XIIa preparation was tested with registered methods and found negative for HIV (1 and 2) antibodies, HBs Ag and HCV antibodies. BSA was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

**USE:**

For *in vitro* use only.

All research studies and protocols where a source of highly purified human Factor XIIa is necessary.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

Restored :

- 24 hours at room temperature (18-25°C).
- 72 hours at 2-8°C.
- 2 months at -20°C or below.

Lyophilized:

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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