Anti-Thrombin (AT) Anti-Thrombine (AT)

Vial of 1.5 mg (PP004B) /Flacon de 1.5 mg (PP004B)
Set of 10 x 1.5 mg (PP004L)/ Coffret de 10 x 1.5 mg (PP004L)

Ref. PP004B / PP004L

STORE AT
2-8°C

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC
PROCEDURES.

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English, last revision 03-2017

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ORIGIN:

Highly purified from human plasma using a combination of affinity and ion-exchange chromatographies.

PRESENTATION:

Vial containing approximately 1.5 mg of AT corresponding to about 10 Plasma Equivalent Units or PEU (1 unit corresponds to the amount of AT in 1 mL human plasma, or about 150µg), which is about 10 International Units (IU) (when determined and standardized against the WHO/NIBSC International standard for Antithrombin, plasma, human).

This kit is for research use only and should not be used for patient diagnosis or treatment.

RECONSTITUTION:

Each vial must be restored with 1 mL distilled water, in order to get a concentration of about 1.5 mg/mL. Dilute in the adequate buffer, if necessary.

EXCIPIENTS:

Content per vial: Glycine (20 mg), Hepes (6 mg), Sodium Chloride (9 mg).

PURITY:

One major band at about 58,000 daltons on SDS-PAGE.

ACTIVITY:

In a clotting assay for the thrombin inhibition capacity in presence of an excess of unfractionated heparin : 10 units AT neutralize at least 2,000 NIH units of thrombin. Tested for anti-Xa activity in presence of heparin, with a chromogenic assay (> 6 IU/mg).

Factor Xa neutralisation capacity is > 500 µg of Xa/mg AT.

The absence of heparin is verified after purification, using a toluidine blue assay (refer to: F.C MacIntosh, "A colorimetric method for the standardization of heparin preparations", Biochem, 1941, 35:776-782), and an anti-Xa specific chromogenic assay.

The exact activity is reported on the analysis certificate for each lot.

VIRAL SAFETY:

The plasma used for AT purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HCV antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

USE:

For *in vitro* use only. All research studies and protocols where a source of highly purified human AT is necessary.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

STABILITY:

Restored :

- at least 72 hours at room temperature (18-25°C).
- at least 1 month at 2-8°C.
- 6 months frozen at -20°C or below.

(tested for anti-Xa chromogenic activity in presence of heparin: > 6 IU/mg).

Lyophilized :

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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