Human Prothrombin

REF PP006B, Pn 1 x 1 mL

REF PP006K, Pn 6 x 1 mL

Human Prothrombin
FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.

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English, last revision 07-2020

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INTENDED USE:

All research studies and protocols where a source of highly purified Prothrombin is necessary.

This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

[Pn] Highly purified prothrombin from citrated human plasma.

Vial containing at least 0.9 mg of highly purified human prothrombin in presence of additives and preservatives.

1 vial of 1 mL (PP006B).

6 vials of 1 mL (PP006K).

- **Excipients:** Glycine, Sodium Chloride, HEPES.
- **Purity:** One major band at 72,000 daltons on SDS-PAGE.
- **Activity:** Tested for its clotting activity in a Factor II deficient plasma system. Prothrombin has about 10 Plasma Equivalent Units (amount in 1 mL of normal human plasma) activity per mg.

WARNING AND PRECAUTIONS:

- Some reagents provided in these kits contain materials of human origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious
- Waste should be disposed of in accordance with applicable local regulations.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- This device of *in vitro* use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Gently remove the freeze-drying stopper, to avoid any product loss when opening the vial.

[Pn] Reconstitute the contents of each vial with exactly **1 mL of distilled water**, shake vigorously until fully dissolved. Homogenize and allow to stabilize for 15 minutes at room temperature (18-25°C), homogenize before use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

[Pn] Reagent stability after reconstitution, free from any contamination or evaporation, and stored closed, is of:

- **7 days** at 2-8°C.
- **78 hours** at room temperature (18-25°C).
- **6 months** frozen at -20°C or less*

*Thaw only once, as rapidly as possible at 37°C, adapting the incubation period to the volume of reagent. The stability of the thawed reagent should be checked under laboratory work conditions.

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

Changes compared to the previous version.

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