# Recombinant Hirudin

## Hirudine recombinante

Vial of  $\approx 16,000$  ATU / flacon de  $\approx 16\ 000$  ATU (RE120A)

Ref: RE120A

STORE AT  
2-8°C

FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC  
PROCEDURES.

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Recombinant protein

### **PRESENTATION:**

Vial containing about 16,000 ATU of rec-Hirudin.

**This specification is based on the hirudin anti-thrombin activity expressed in ATU\*.**

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

### **RECONSTITUTION:**

Restore the vial with 1 ml of distilled water. When restored, the solution contains a rec-Hirudin solution with an activity of approximately 16,000 ATU\*/ml.

The measured recombinant hirudin activity is given for each lot on the Certificated of Analysis.

### **EXCIPIENTS:**

- Glycine, Hepes, Sodium Chloride, **pH 7.50**

### **ACTIVITY:**

Hirudin has an activity of about 16,000 ATU\*/vial

### **USE:**

For *in vitro* use only.

For research use only (RUO)

Thrombin inhibitor, which can be used in any *in vitro* laboratory research application where an efficient inhibition of Thrombin is required.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

### **STABILITY:**

Restored:

- 2 days at room temperature
- 7 days at 2-8°C
- Frozen at -20°C or below\*.

\*Thaw once as rapidly as possible at 37°C, adapt duration to the volume of reagent. The stability of the thawed reagent should be verified in the working conditions of the user laboratory.

Lyophilized :

- Until the expiration date printed on the vial.

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