PLASMA PF4 CONTROL
SC013KFOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSIS PROCEDURES.POUR LA RECHERCHE UNIQUEMENT.
NE PAS UTILISER DANS LES PROCEDURES DE DIAGNOSTIC.**Intended use**

Lyophilised plasmas, at a defined Platelet Factor 4 (PF4) concentration, as measured with the ZYMUTEST-PF4 kit (ref RK006A), for the quality control of PF4 antigenic assays. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

Reagents:

Each kit contains:

- 3 vials of 0.5 ml of normal plasma, lyophilised, at a high PF4 concentration (Control I)
- 3 vials of 0.5 ml of normal plasma, lyophilised, at a low PF4 concentration (Control II).

The PF4 concentrations and the acceptance ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

Reagent Preparation and Stability:

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

Excipients:

Glycine, stabilizers.

Preparation:

- **PF4 Control I (high):** restore the vial with 0.5 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **PF4 Control II (low):** restore the vial with 0.5 ml distilled water ; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These plasma are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota : Plasma used for PF4 controls were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Application

The PF4 concentration, indicated on the flyer, for Control I and Control II, has been measured with the ZYMUTEST PF4 (RK006A), Elisa kit.

This assay is standardised by reference to the NIBSC international standard.

The target values and the acceptance ranges, indicated on the flyer, are valid for ZYMUTEST PF4 only. The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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