PLASMA PROTEIN S CONTROL
SC050K**FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSIS PROCEDURES.****POUR LA RECHERCHE UNIQUEMENT.
NE PAS UTILISER DANS LES PROCEDURES DE DIAGNOSTIC.****Intended use**

Lyophilised plasmas, at a defined Protein S concentration, as measured either with the ZYMUTEST-Free Protein S kit (ref RK015A) or the ZYMUTEST-Total Protein S kit (ref RK021A), for the quality control of Protein S antigenic assays. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

Reagents:

Each kit contains:

- 3 vials of 0.5 ml of normal plasma, lyophilised, at a high Protein S concentration (Control I)
- 3 vials of 0.5 ml of normal plasma, lyophilised, at a low Protein S concentration (Control II).

The Protein S concentrations and the acceptance ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

Reagent Preparation and Stability:

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

Excipients:

Glycine, stabilizers.

Preparation:

- **Protein S Control I (high):** restore the vial with 0.5 ml distilled water; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **Protein S Control II (low):** restore the vial with 0.5 ml distilled water; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These plasmas are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota : Plasmas used for Protein S controls were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Application

The Protein S concentrations, indicated on the flyer, for Control I and Control II, have been measured with the ZYMUTEST Free Protein S (ref RK015A) or ZYMUTEST Total Protein S (ref RK021A), Elisa kits.

These assays are standardised by reference to the NIBSC international standard.

The target values and the acceptance ranges, indicated on the flyer, are valid for ZYMUTEST Protein S assays (Free and Total). The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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