# BIOPHEN CS-01(38)

## Thrombin Chromogenic substrate

Vial of 25 mg/flacon de 25 mg  
Set of 12x25 mg (229001C) / Coffret de 12x25 mg (229001C)  
Set of 50x25 mg (229001Z) / Coffret de 50x25 mg (229001Z)

Ref. 229001 / 229001C / 229001Z

**STORE AT**  
**2-8 °C**

**FOR RESEARCH USE ONLY.**  
**NOT FOR USE IN DIAGNOSTIC**  
**PROCEDURES.**

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English, last revision 03-2017

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**ORIGIN:**

Synthetic generic Chromogenic substrate, highly purified and stabilized.

**SPECIFICITY:**

Recommended substrate for Thrombin.

**PRESENTATION:**

Each vial contains about 25 mg of BIOPHEN CS-01(38), lyophilised in presence of Mannitol as a bulking agent.

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**RECONSTITUTION:**

According to the research protocol used, the BIOPHEN CS-01(38) chromogenic substrate can be restored with variable volumes of distilled water; for example 5 mL can be used for a substrate concentration of 5 mg/mL, or 20 mL for a substrate concentration of 1.25 mg/mL.

**PEPTIDE SEQUENCE:** H-D-Phe-Pip-Arg-pNa, 2HCl.

**CHEMICAL STRUCTURE:** C<sub>27</sub>H<sub>36</sub>N<sub>8</sub>O<sub>5</sub>, 2HCl.

**MOLARITY:** 45 µMoles.

**PURITY GRADE:** ≥ 95%.

**MOLECULAR WEIGHT (-2HCL):** 553 Da.

**MOLECULAR WEIGHT (+2HCL):** 625 Da.

**FREE PNA CONTENT (MEASURED ON SUBSTRATE AT 2.5G/L):** OD<sub>405nm</sub> ≤ 0.30.

**USE:**

All research studies and protocols where a source of chromogenic substrate for Thrombin is required.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**RESPECTIVE REACTIVITIES:**

Thrombin	FXa	Plasmin	Kallicrein	aPC
100	5	5	60	40

Assay conditions must be established for making the substrate totally specific for Thrombin.

**STABILITY:**

Restored:

- 7 days at room temperature (18-25°C).
- 3 months at 2-8°C.
- Do not freeze.

Lyophilized:

- Until the expiration date printed on the vial.

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