BIOPHEN CS-21(66)

Activated Protein C Chromogenic substrate

Vial of 25 mg/flacon de 25 mg/Provetta da 25 mg

Ref. 229021

RUO

ENGLISH

STORE AT 2-8 °C For Research use only
Not for use in diagnosis procedures.

Origin: Synthetic generic Chromogenic substrate, highly purified and stabilized.

Specificity: Recommended substrate for Activated Protein C.

Presentation: Vial containing 25 mg of BIOPHEN CS-21(66), lyophilised in presence of Mannitol as a bulking agent.

Reconstitution: According to the research protocol used, the BIOPHEN CS-21(66) chromogenic substrate can be restored with variable volumes of distilled water; for example 5 mL can be used for a substrate concentration of 5 mg/mL, or 20 mL for a substrate concentration of 1.25 mg/mL.

Peptide sequence: p-Glu-Pro-Arg-pNa. HCl

Chemical structure: C₂₂H₃₀N₆O₆, HCl

Molarity: ≈50 µmol/vial

Purity grade: > 95%.

Molecular weight: 502.5.

Free pNa content (measured on substrate at 2.5g/L): OD_{405nm} ≤ 0.30

Use:

All research studies and protocols where a source of chromogenic substrate for Activated Protein C is required. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

Respective reactivities:

APC	FXa	Plasmin	Kallicrein	Thrombin
100	7	85	150	125

Assay conditions must be established for making the substrate totally specific for activated protein C.

Stability:

Restored :

- 7 days at room temperature
- 3 months at 2-8 °C
- **Do not freeze.**

Lyophilized :

Until the expiration date printed on the vial.

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