Tris-NaCl buffer, with Albumin – pH 7.40

Tampon Tris-NaCl Albuminé – pH 7,40

Vial of 50 mL (AR005A) / Flacon de 50 mL (AR005A)
Set of 3x50 mL (AR005K) / Coffret de 3x50 mL (AR005K)

Ref. AR005A / AR005K

STORE AT
2-8°C**FOR RESEARCH USE ONLY.**
NOT FOR USE IN DIAGNOSTIC
PROCEDURES**UNIQUEMENT A USAGE DE**
RECHERCHE. NE PAS UTILISER A DES
FINS DE DIAGNOSTICS.**CONSERVER A**
2-8°C**Lot:****Expiration:**

English, last revision 06-2017

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COMPOSITION:

Buffer containing Tris (0.05M) (Tris (hydroxymethyl) aminomethane), Sodium Chloride (0.15M), Bovine Serum Albumin (BSA at 1%) and 0.90 g/L Sodium Azide (as preservative).

*CAUTION: Sodium azide (NaN₃), may react with lead and copper plumbing to form highly explosive metal azides. Flush with large volumes of water when discarding into a sink.

This kit is for research use only and should not be used for patient diagnosis or treatment.

PH:

7.40.

PRESENTATION:

Vial of 50 mL, ready to use.

USE:For *in vitro* use only.

Dilution buffer for biological tests, recommended mainly for diluting purified proteins and enzymes (Albumine prevents adsorption of purified proteins on recipient walls).

Can be used as "buffer A" recommended by the European Pharmacopoeia for Heparin anti-IIa assay.

Note 1: Take care of stabilizing the solution at room temperature (18-25°C) for 30 minutes before use.

Note 2: Refer to the specific associated kit inserts for instructions concerning the tests to be performed.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

STABILITY:

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

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