FPA anticoagulant solution

Solution anticoagulante pour FPA

(Vial of 20 mL / Flacon de 20 mL)

Ref. AR013A

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English, last revision 06-2017

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COMPOSITION:

Solution containing TriSodium citrate, Heparin, Aprotinin, Hirudin, and Sodium azide (0.9 g/L) as preservative.

***CAUTION: Sodium azide (NaN₃), may react with lead and copper plumbing to form highly explosive metal azides. In order to avoid this risk, flush with large volumes of water when discarding into a sink.**

This kit is for research use only and should not be used for patient diagnosis or treatment.

PRESENTATION:

20 mL vial, ready to use.

USE:

For *in vitro* use only.

Special anticoagulant solution (ACS) for the assay of FPA.

Ready to use

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

SPECIMEN COLLECTION AND PREPARATION:

Blood (9 vol.) must be collected by a net venipuncture on the special anticoagulant mixture (1 vol.) provided for the assay of FPA (contains Trisodium Citrate, Heparin, Hirudin, Aprotinin and sodium azide) and discarding the first drops. Blood must be rapidly mixed with the anticoagulant and centrifuged at 2,500 g for 20 min. Plasma supernatant must be decanted and is then ready for bentonite treatment. Plasma must be treated with bentonite within 8 hours after collection or deep frozen within 4 hours, and stored up to 1 month at -20°C or below. Just before use, it must be thawed for 30 min. at 37°C in a water-bath, and treated with bentonite.

STABILITY:

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

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