Para-Nitro-Aniline (pNA)

(Vial of 5 ml / Flacon de 5 ml)

Ref. AR014A

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Lot:

Expiration:

Product:

p-Nitroanilin solution (pNA) .

Standard solution, with an absorbance at 405 nm, for controlling photometers and spectrophotometers used for chromogenic assays, based on release of pNA.

Composition:

pNA solution at 0.10 mmol/L.

Use:

For in vitro use only.

The pNA standard solution allows controlling the accuracy of spectrophotometers and photometers (at 405 nm) used for measuring the yellow colour developed in chromogenic assays where pNA is released.

The target value for absorbance is:

- Using 1 cm cuvette:
A405 = 0.965 ± 0.015

INDICATIVE: For Micro ELISA wells with 200 µl reagent :
A405 = 0.530 ± 0.030

Protocol:

Photometer cuvette, 1 cm optical path	Microwell (INDICATIVE)
Set up the spectrophotometer at 405 nm, and let it to equilibrate	Fill the wells either with 200 µl of distilled water (1 strip of 8 wells) or pNA standard solution (another strip of 8 wells).
Introduce 800 µl of pNA standard in a cuvette.	
Introduce 800 µl of distilled water.	
Measure the absorbance at 405 nm.	Measure the mean value A405, and subtract the mean A405 value obtained for distilled water from the one measured for the standard pNA solution
Subtract the absorbance of water from that of pNA for obtaining the net absorbance of the pNA solution.	

Results:

If the measured absorbance is out of the acceptance range, correct the A405 measured for the chromogenic assays by a correcting factor obtained by dividing the target A405 value of the pNA standard solution by the measured A405 value.

Stability and storage of the vial:

- Until the expiration date printed on the label, when the vial is stored at 2-8 °C
- Store the vial protected from bright sunlight.

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