# Fibrinolysis Sample Diluent (for ELISA kit)

(Vial of 50 mL / Flacon de 50 mL)

Ref. AR016A

**STORE AT**  
**2-8 °C**

**FOR RESEARCH USE ONLY.**  
**NOT FOR USE IN DIAGNOSTIC**  
**PROCEDURES.**

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**Expiration:**

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**COMPOSITION:**

Phosphate buffer saline (9 g/L NaCl) containing bovine serum albumin (1%) and Tween 20 (0.1%).  
Contains Kathon CG as preservative.  
Reagent ready to use.

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**VIRAL SAFETY:**

BSA was prepared from bovine plasma; which was tested for the absence of infectious agents, and collected from animals free from BSE. However, no test may totally exclude the absence of infectious agents. As any product of bovine origin, this reagent must be used with all the cautions required for handling a material potentially infectious.

**PH:**

7.50.

**USE:**

For *in vitro* use only.  
Dilution buffer for fibrinolysis investigations by ELISA method.  
Ready to use.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.  
When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

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