# Autoimmunity Sample Diluent (for ELISA kit)

(Vial of 50 mL / Flacon de 50 mL)

Ref. AR018A

**STORE AT**  
**2-8 °C**

**FOR RESEARCH USE ONLY.**  
**NOT FOR USE IN DIAGNOSTIC**  
**PROCEDURES.**

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**Expiration:**

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**COMPOSITION:**

Phosphate buffer saline (9 g/L NaCl) containing goat serum (<20%) and Tween 20 (0.1%).

Contains Sodium Azide (NaN<sub>3</sub>) as preservative.

Reagent ready to use.

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**VIRAL SAFETY:**

Goat serum: No test may totally exclude the absence of infectious agents. As any product of biological origin, this reagent must be used with all the cautions required for handling a material potentially infectious.

**WARNING:**

The Autoimmunity Sample Diluent contains Sodium Azide (NaN<sub>3</sub>), which may react with lead and copper plumbing to form highly explosive metal azides. Flush with large volumes of water when discarding into a sink.

**PH:**

7.50.

**USE:**

For *in vitro* use only.

Dilution buffer for testing autoantibodies by ELISA method.

Ready to use.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

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