Purified bovine Fibrinogen

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Vial of 100 mg/flacon de 100 mg

Ref. BP101A

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C**Origin:**

Highly purified from citrated bovine plasma

Presentation:

Vial containing approximately 100 mg of highly purified bovine fibrinogen.

Reconstitution:

Each vial must be restored with 5 mL distilled water. The solution can be diluted to the desired concentration when required.

Excipients:

- Glycine, trisodium citrate and HEPES buffer
- 45 mg sodium chloride/vial

Purity:One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clottability $\geq 98\%$.**Activity:**Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/mL (100 μ L), is clotted by bovine thrombin, at 5 NIH/mL (100 μ L), the clotting time must be ≤ 12 sec.**Viral safety:**

The plasma used for fibrin purification was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

Use:For *in vitro* use only.
All research studies and protocols where a source of highly purified bovine fibrinogen is necessary.**Stability:**

- Restored :**
- 48 hours at room temperature (18-25 °C)
 - 48 hours at 2-8 °C

Lyophilized : Until the expiration date printed on the vial.

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CONSERVER A
2-8 °C

ANALYSIS CERTIFICATE

ANALYTICAL DATA	SPECIFICATIONS
1. Protein Content (Lowry method; A280 nm) mg/vial	> 90 mg/vial
2. SDS-PAGE (5 % acrylamide) 1 single band of about daltons	1 major band of about 340,000 daltons
3. Clotting time (5 NIH/mL Thrombin) sec	≤12 sec
4. Coagulability %	≥ 98%

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DATE :

CONCLUSIONS : Passed/Conforme Refused/Non-conforme

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