Purified bovine Factor X

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Vial of 100 µg/Flacon de 100 µg
Set of 6x100 µg/Coffret de 6x100 µg

Ref. BP102A / BP102K

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English, last revision 12-2016

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ORIGIN:

Highly purified from citrated bovine plasma.

PRESENTATION:BP 102A: 1 vial containing approximately 100 µg of highly purified bovine Factor X.
BP102K: Set of 6 vials containing approximately 100 µg of highly purified bovine Factor X.**This kit is for research use only and should not be used for patient diagnosis or treatment.****RECONSTITUTION:**

Each vial must be restored with 1 mL distilled water. The solution can be diluted to the desired concentration, if required.

EXCIPIENTS:Glycine and Hepes.
9 mg sodium chloride/vial.**PURITY:**

One major band of 55,000 daltons on SDS-PAGE.

ACTIVITY:

Tested for :

- a. Clotting activity (Factor X deficient plasma) : Activity is about 80 U/mg
b. Amidolytic activity, following activation with RVV.
(1 Unit "U" is the amount in 1 mL normal bovine plasma).

VIRAL SAFETY:

The plasma used for Factor X purification was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

USE:

For in vitro use only.

All research studies or protocols where a source of highly purified bovine Factor X is necessary.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**STABILITY:**

Restored : - 24 hours at room temperature (18-25°C).
- 72 hours at 2-8 °C.
- 6 months at -20°C or below.
Lyophilized : Until the expiration date printed on the kit.

The stability of the reconstituted reagent should be checked under laboratory work conditions.

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