Bovine Plasminogen

Plasminogène bovin

Vial of 5 mg/flacon de 5 mg

Ref. BP105C

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English, last revision 03-2017

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ORIGIN:

Highly purified from bovine plasma using affinity chromatography.

PRESENTATION:

Vial containing approximately 5 mg of plasminogen corresponding to about 20 Plasma Equivalent Units or PEU (amount in 1 mL of normal bovine plasma).

This kit is for research use only and should not be used for patient diagnosis or treatment.**RECONSTITUTION:**

Each vial must be restored with 5 mL distilled water, in order to obtain a solution at 1 mg/mL.

EXCIPIENTS:

- Phosphate buffer at pH 7.50
- 4.5 mg sodium chloride/vial.

PURITY:

One major band on SDS-PAGE.

SPECIFIC ACTIVITY:

About 4 Plasma Equivalent Unit/mg as tested with a chromogenic assay following activation to plasmin by urokinase, in presence of EACA.

VIRAL SAFETY:

The plasma used for plasminogen purification was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

USE:For *in vitro* use only.

All research studies and protocols where a source of highly purified bovine plasminogen is necessary.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**STABILITY:**

Restored:

- 24 hours at room temperature (18-25°C).
- 72 hours at 2-8 °C.
- 2 months frozen at -30°C or below.

Lyophilized:

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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