# Urokinase type Plasminogen Activator (uPA)

Vial of 1,000 units/flacon de 1 000 unités

Ref. EZ005A

**STORE AT**  
**2-8 °C**

**FOR RESEARCH USE ONLY.**  
**NOT FOR USE IN DIAGNOSTIC**  
**PROCEDURES.**

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English, last revision 03-2017

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**ORIGIN:**

Purified Urokinase, extracted from human urine.  
The preparation contains about 85% High Molecular Weight Urokinase and 15% Low Molecular Weight Urokinase.

**PRESENTATION:**

Vial containing approximately about 1,000 Units of Urokinase lyophilised in presence of additives and preservatives.  
**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**RECONSTITUTION:**

Each vial must be restored with 1 mL distilled water in order to yield a preparation containing about 1,000 IU/mL.

**EXCIPIENTS:**

Sodium Phosphate, Sodium Chloride, Poly-Ethylene-Glycol 6000 (PEG 6000), Bovine serum Albumine and Mannitol.

**PURITY:**

Purified Urokinase with a specific activity  $\geq 75,000$  IU/mg of protein.  
One major band of 54,000 daltons corresponding to High Molecular Weight Urokinase and one minor one of 33,000 daltons corresponding to Low Molecular Weight Urokinase, on SDS-PAGE.

**VIRAL SAFETY:**

The Urokinase preparation contains BSA, which was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

**USE:**

For *in vitro* research use only.  
Used in any application where a source of Purified Urokinase is required.  
Used for the assay of human plasminogen.  
Activates plasminogen in order to generate Plasmin.  
**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

Restored :  
- 24 hours at room temperature (18-25°C).  
- 72 hours at 2-8°C.  
- 2 months at -20°C or below.  
Lyophilized :  
- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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DATE :

CONCLUSIONS :  Passed/*Conforme*  Refused/*Non-conforme*

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