Urokinase type Plasminogen Activator (uPA)

Vial of 100,000 units/flacon de 100 000 unités

Ref. EZ005C

STORE AT
2-8 °C

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC
PROCEDURES.

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English, last revision 03-2017

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ORIGIN:

Purified Urokinase, extracted from human urine.
The preparation contains about 85% High Molecular Weight Urokinase and 15% Low Molecular Weight Urokinase.

PRESENTATION:

Vial containing approximately about 100,000 Units of Urokinase lyophilized in presence of additives and preservatives.

This kit is for research use only and should not be used for patient diagnosis or treatment.

RECONSTITUTION:

Each vial must be restored with 0.5 mL distilled water in order to yield a preparation containing about 200,000 IU/mL (100,000 IU/vial).

EXCIPIENTS:

Sodium Phosphate, Sodium Chloride.
Excipients are limited to the strict necessary in order to allow using the Urokinase preparation in the required buffer. Following reconstitution with 0.5 mL distilled water, a solution containing 200,000 IU/mL is obtained. The required buffer can be added directly into the vial in order to obtaining the stock or working solution.

PURITY:

Purified Urokinase with a specific activity $\geq 75,000$ IU/mg of protein.
One major band of 54,000 daltons corresponding to High Molecular Weight Urokinase and one minor one of 33,000 daltons corresponding to Low Molecular Weight Urokinase, on SDS-PAGE.

VIRAL SAFETY:

Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

USE:

For *in vitro* research use only.
Used in any application where a source of Purified Urokinase is required.
Used for the assay of human plasminogen.
Activates plasminogen in order to generate Plasmin.
The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

STABILITY:

Restored :

- 24 hours at room temperature (18-25°C).
- 72 hours at 2-8°C.
- 2 months at -20°C or below.

Lyophilized :

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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