Streptokinase

Vial of 100,000 units/flacon de 100 000 unités

Ref. EZ008C

STORE AT
2-8 °C

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC
PROCEDURES.

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Purified Streptokinase.

PRESENTATION:

Vial containing approximately about 100,000 Units of Streptokinase lyophilised in presence of additives and preservatives.

This kit is for research use only and should not be used for patient diagnosis or treatment.

RECONSTITUTION:

Each vial must be restored with 1 mL distilled water.

EXCIPIENTS:

Sodium Phosphate, Sodium Chloride, Poly-Ethylene-Glycol 6000 (PEG 6000), Bovine serum Albumine.

PURITY:

Purified streptokinase with a specific activity between 15,000 to 25 000 IU/mg of protein.

VIRAL SAFETY:

The streptokinase preparation contains BSA, which was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

USE:

For *in vitro* use only.

Forms complexes with human plasminogen in order to generate a Plasmin activity. Used for the assay of human plasminogen.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

STABILITY:

Restored :

- 24 hours at room temperature (18-25°C).
- 72 hours at 2-8°C.
- 2 months at -20°C or below.

Lyophilized :

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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