# Purified human Factor IXa

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Vial of 100 IU / flacon de 100 UI

Ref. EZ010B

**STORE AT**  
**2-8 °C**

For research use only

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**CONSERVER A**  
**2-8 °C****Origin:**

Highly purified Factor IXa prepared from Factor IX extracted from a citrated human plasma pool, activated with agarose bound Factor XIa. MW of about 45,000 Da.

**Presentation:**

Vial containing approximately 100 IU of highly purified human factor IXa in presence of additives and preservatives.  
(1 International Unit of FIXa corresponds to the Factor IXa determined against the International Standard for Factor IXa, WHO/NIBSC).

**Reconstitution:**

Each vial must be restored with 1 ml distilled water, in order to get a concentration of about 100 IU/ml. It can be diluted to the desired concentration, when required, in the adequate buffer.

**Excipients:**

Content per vial: Hepes (6mg), Mannitol (2%), Prionex (traces), Sodium chloride (9mg), at pH 7.50.

**Purity:**

The Factor IX zymogen used for preparing Factor IXa has one major band of 55,000 daltons on SDS-PAGE.

**Activity:**

Tested for its amidolytic activity using the BIOPHEN Factor IXa kit (Ref 221812) by reference to the NIBSC/WHO International Standard in force.

**Viral safety:**

The human plasma used for Factor IXa preparation was tested with registered methods and found negative for HIV (1 and 2) antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For *in vitro* use only. All research studies and protocols where a source of highly purified human Factor IXa is necessary. Studies on heparin and Low Molecular weight heparin.

**Stability:**

**Restored :** - 24 hours at room temperature (18-25 °C)  
- 48 hours at 2-8 °C  
- 6 months at -20 °C or below.

**Lyophilized :** Until the expiration date printed on the vial.

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