# Purified Human Fibrinogen

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Vial of 5 mg/flacon de 5 mg

Ref. PP001A

**STORE AT**  
**2-8 °C**

For research use only

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Highly purified from citrated human plasma

**Presentation:**

Vial containing approximately 5 mg of highly purified human fibrinogen.

**Reconstitution:**

Each vial must be restored with 1 mL distilled water. It can be diluted to the desired concentration when required.

**Excipients:**

- Glycine, trisodium citrate and HEPES buffer
- 6 mg sodium chloride/vial

**Purity:**One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clottability  $\geq 98\%$ .**Activity:**Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/mL (200  $\mu$ l), is clotted by thrombin at 5 NIH/mL (100  $\mu$ l) the clotting time must be below 15 sec.**Viral safety:**

The plasma used for fibrinogen purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**For *in vitro* use only.  
All research studies and protocols where a source of highly purified fibrinogen is necessary.**Stability:**

- Restored :**
- 48 hours at room temperature (18-25 °C)
  - 96 hours at 2-8 °C

**Lyophilized :** Until the expiration date printed on the vial.

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CONSERVER A  
2-8 °C

### ANALYSIS CERTIFICATE

| ANALYTICAL DATA                                                                                                | SPECIFICATIONS                           |
|----------------------------------------------------------------------------------------------------------------|------------------------------------------|
| <b>1. Protein Content (Lowry method)</b><br><p style="text-align: center;">..... mg/vial</p>                   | > 4.5 mg/vial                            |
| <b>2. SDS-PAGE (5 % acrylamide)</b><br><p style="text-align: center;">1 single band of about ..... daltons</p> | 1 major band of about<br>340,000 daltons |
| <b>3. Clotting time (5 NIH/mL Thrombin)</b><br><p style="text-align: center;">..... sec</p>                    | < 15 sec                                 |
| <b>4. Coagulability</b><br><p style="text-align: center;">..... %</p>                                          | ≥ 98%                                    |

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DATE :

CONCLUSIONS :  Passed/Conforme  Refused/Non-conforme

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