**Purified Human Fibrinogen, Salt Free**  
**Fibrinogène humain purifié sans sels**  
Vial of 100 mg/flacon de 100 mg

Ref. PP001R

**STORE AT**  
**2-8 °C**

For research use only

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**CONSERVER A**  
**2-8 °C****Origin:**

Highly purified fibrinogen extracted from citrated human plasma.

**Presentation:**

Vial containing approximately 100 mg of highly purified human fibrinogen, salt free.

**Reconstitution:**

Each vial must be restored with 10 mL of the desired working buffer in order to obtain a 10 mg/mL of fibrinogen solution. Incubate preferentially at 37 °C, until complete dissolution. Shake gently, avoiding formation of foam. Total dissolution of fibrinogen content may require 1 hour or more. This solution can be diluted to the desired concentration when required.

**Excipients:**

This fibrinogen preparation is salt free and excipient free.

**Purity:**

One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clottability ≥ 98%.

**Activity:**

Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/mL (100 µL), is clotted by thrombin, at 5 NIH/mL (100 µL), the clotting time must be below 15 sec.

**Viral safety:**

The plasma used for fibrinogen purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For in vitro use only.  
All research studies and protocols where a source of highly purified fibrinogen is necessary.

**Stability:**

**Restored :** - 48 hours at room temperature (18-25 °C)  
- 96 hours at 2-8 °C

**Lyophilized :** Until the expiration date printed on the vial.

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