# Human platelet factor 4 (PF4)

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Vial of 100 µg/flacon de 100 µg

Ref. PP003A

**STORE AT**  
**2-8 °C**

For research use only

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**CONSERVER A**  
**2-8 °C****Origin:**

Human platelet concentrates.

**Presentation:**

Vial containing approximately 100 µg of PF4 in the tetrameric form.

**Reconstitution:**

Each vial must be restored with 1 mL distilled water.  
When restored, the solution contains about 100 µg/mL of PF4 in a 0.5 M sodium chloride buffer.

**Nota** : A high salt concentration ( $\geq 0.5$  M) is recommended for long term storage ( $> 6$ h) after reconstitution.

**Excipients:**

– Glycine and HEPESbuffer - 30 mg sodium chloride/vial

**Purity:**

One major band of 8,000 daltons on SDS-PAGE, for reduced PF4 monomer.

**Activity:**

Tested in a clotting assay for the heparin inhibition: 1 vial (100 µg) neutralizes at least 1.25 IU of heparin.

**Viral safety:**

The material used for PF4 preparation was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For in vitro use only.  
All research studies and protocols where a source of highly purified PF4 is necessary.

**Stability:**

**Restored** : - 8 hours at room temperature (18-25 °C)  
- 24 hours at 2-8 °C  
- 2 months frozen at -30 °C or below

**Lyophilized** : Until the expiration date printed on the vial.

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