# Anti-Thrombin (AT)

# Anti-Thrombine (AT)

Vial of 10 mg/flacon de 10 mg

Ref. PP004C

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English, last revision 03-2017

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**ORIGIN:**

Purified from human plasma using a combination of affinity chromatographies.

**PRESENTATION:**

Vial containing approximately 10 mg of AT corresponding to about 60 Plasma Equivalent Units or PEU (1 unit corresponds to the amount of AT in 1 mL human plasma, or about 150 µg), which is about 60 International Units (IU) (when determined and standardized against the WHO/NIBSC International standard for Antithrombin, plasma, human).

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**RECONSTITUTION:**

Each vial must be restored with 2 mL distilled water, in order to get a concentration of about 5 mg/mL, in a physiological phosphate saline buffer (pH 7.5). Dilute in the adequate buffer, if necessary.

**EXCIPIENTS:**

Content per vial: Phosphate (16 mg), Sodium Chloride (18 mg).

**PURITY:**

One major band at about 58,000 daltons on SDS-PAGE.

**ACTIVITY:**

In a clotting assay for the thrombin inhibition capacity in presence of an excess of unfractionated heparin: 10 mg AT neutralize at least 12,000 NIH units of thrombin. Tested for anti-Xa activity in presence of heparin, with a chromogenic assay (> 6 IU/mg).

Factor Xa neutralisation capacity is > 500 µg of Xa/mg AT.

The absence of heparin is verified after purification, using a toluidine blue assay (refer to: F.C MacIntosh, "A colorimetric method for the standardization of heparin preparations", Biochem, 1941, 35:776-782), and an anti-Xa specific chromogenic assay.

The exact activity is reported on the analysis certificate for each lot.

**VIRAL SAFETY:**

The plasma used for AT purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HCV antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**USE:**

For *in vitro* use only. All research studies and protocols where a source of highly purified human AT is necessary.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**STABILITY:**

Restored :

- 72 hours at room temperature (18-25°C).
- 7 days at 2-8°C.
- 2 months frozen at -30°C or below.

Lyophilized :

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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