# Plasminogen Plasminogène

Vial of 1.5 mg/flacon de 1.5 mg

Ref. PP005B

**STORE AT**  
**2-8 °C**

For research use only

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**CONSERVER A**  
**2-8 °C****Origin:**

Highly purified from human plasma using affinity chromatography.

**Presentation:**

Vial containing approximately 1.5 mg of plasminogen corresponding to about 6 Plasma Equivalent Unit or PEU (amount in 1 mL of normal human plasma).

**Reconstitution:**

Each vial must be restored with 2 mL distilled water, in order to obtain a solution at 0.75 mg/mL.

**Excipients:**

- Phosphate buffer at pH 7.50
- 9 mg sodium chloride/vial, Mannitol

**Purity:**

One major band at 88,000 daltons on SDS-PAGE

**Specific activity:**

About 4 Plasma Equivalent Unit/mg as tested with a chromogenic assay following complex formation with streptokinase.

**Viral safety:**

The plasma used for plasminogen purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For *in vitro* use only.

All research studies and protocols where a source of highly purified plasminogen is necessary.

**Stability:**

- Restored :**
- 24 hours at room temperature (18-25 °C)
  - 72 hours at 2-8 °C
  - 6 months frozen at -30 °C or below

**Lyophilized :** Until the expiration date printed on the vial.

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