# Human Prothrombin

## Prothrombine humaine

Vial of 1 mg (PP006B) /flacon de 1 mg (PP006B)  
Set of 6x1 mg (PP006K) /Coffret de 6x1 mg (PP006K)

Ref. PP006B / PP006K

**STORE AT**  
**2-8 °C**

For *in vitro* research use only

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**CONSERVER A**  
**2-8 °C**

**Origin:**

Highly purified prothrombin from citrated human plasma.

**Presentation:**

Vial containing approximately 1 mg of highly purified human prothrombin in presence of additives and preservatives.

**Reconstitution:**

Each vial must be restored with 1 ml distilled water.

**Excipients:**

– Glycine, Sodium Chloride, HEPES.

**Purity:**

One major band at 72,000 daltons on SDS-PAGE

**Specific activity:**

Tested for its clotting activity in a factor II deficient plasma system. Prothrombin has about 10 Plasma Equivalent Units (amount in 1 ml of normal human plasma) activity per mg.

**Viral safety:**

The plasma used for prothrombin purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For *in vitro* use only.

All research studies and protocols where a source of highly purified Prothrombin is necessary.

**Stability:**

**Restored :**

- 78 hours at room temperature (18-25 °C)
- 7 days at 2-8 °C
- 6 months frozen at -20 °C or below

**Lyophilized :** Until the expiration date printed on the vial.

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