Human Protein C

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Vial of 100 µg/flacon de 100 µg

Ref. PP009A

STORE AT
2-8°C

For research use only

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CONSERVER A
2-8°C**Origin:**

Highly purified Protein C from citrated human plasma.

Presentation:

Vial containing approximately 100 µg of highly purified human Protein C in presence of additives and preservatives.

Reconstitution:

Each vial must be restored with 1 ml distilled water.

Excipients:

– Glycine, sodium chloride, HEPES.

Purity:

One major band at 62,000 daltons on SDS-PAGE

Specific activity:

Tested on SDS-PAGE, for its anticoagulant activity, for the antigen concentration and for the absence of factor II in a factor II deficient plasma system.

≥ 20 PEU/vial (1 Plasma Equivalent Unit or PEU is the amount in 1 ml of normal human plasma)

Viral safety:

The plasma used for Protein C purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:For *in vitro* use only.

All research studies and protocols where a source of highly purified Protein C is necessary.

Stability:**Restored :**
- 24 hours at room temperature (18-25°C)
- 3 days at 2-8 °C
- 6 months frozen at -20°C or below**Lyophilized :** Until the expiration date printed on the vial.**Origine :**

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