Prekallikrein Pool

Pool de Prékallikréine

Vial of 2 mL/flacon de 2 mL or Set of 6x2mL/Coffret de 6x2mL

Ref. PP501B / PP501K

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C**Origin:**

Prekallikrein pool from human plasma according to the European Pharmacopeia.

Presentation:

Each Vial contains 2 mL of Prekallikrein pool extracted from human citrated plasma.

Reconstitution:

Each vial must be restored with 2 mL distilled water in order to obtain a prekallikrein pool according to the European Pharmacopeia. Incubate until complete dissolution. Shake gently, avoiding formation of foam. This solution can be diluted to the desired concentration if required.

Excipients:

- About 17 mg sodium chloride/vial
- Mannitol as bulking agent, Tris, sodium chloride

Activity:

Tested with a chromogenic assay using a specific kallikrein substrate, following activation with PK activator (mainly Factor XIIa), established by reference to the NIBSC International Standard.

Viral safety:

The plasma used was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.
All research studies and protocols where a source of Prekallikrein pool is necessary.

Stability:

- Restored :**
- 24 hours at room temperature (18-25 °C)
 - 48 hours at 2-8 °C
 - 2 months frozen at -25 °C or below

Lyophilized : Until the expiration date printed on the vial.

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