# Recombinant human PAI-1 (native)

## PAI-1 humain recombinant (natif)

Vial of 10 µg/flacon de 10 µg

Ref. RP006A

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English, last revision 03-2017

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**ORIGIN:**

Recombinant human PAI-1 (wild type) expressed in E. Coli.

**PRESENTATION:**

Vial containing 10 µg of recombinant human PAI-1 (expressed in E. Coli), lyophilized.

**This kit is for research use only and should not be used for patient diagnosis or treatment.****RECONSTITUTION:**

Restore the vial with 0.2 mL of distilled water, then add 0.8 mL of the appropriate buffer, in order to get a PAI-1 solution at 10 µg/mL.

**Nota :** The dilution buffer must contain Bovine Serum Albumin (i.e. 1%) or any other carrier protein or polymer in order to avoid loss of the protein by adsorption on the container walls.**EXCIPIENTS:**Phosphate, Arginine, Sodium Chloride (0.15 M), pH 7.50.  
Bovine Serum Albumin 1%.**VIRAL SAFETY:**

The bovine plasma used for BSA preparation was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.

**ACTIVITY:**

Recombinant human PAI-1 has an activity of about 50% (tested for its binding to tPA). PAI-1 activity is labile. Recombinant human PAI-1 must be used rapidly or stored frozen in small aliquots and thawed just before use.

**CHARACTERISTICS:**PAI-1 has a molecular weight of about 50,000 daltons.  
Recombinant human PAI-1 preparation has a purity grade > 95%.**USE:**For *in vitro* use only.

All research studies and protocols where a source of recombinant human PAI-1 is necessary.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.****STABILITY:**

Restored :

- 4 hours at room temperature (18-25°C).
- 24 hours at 2-8°C.
- 6 months at -20°C or below.

Lyophilized :

- Until the expiration date printed on the vial.

Freezing conditions and stability of the thawed product should be checked in the working conditions of the laboratory user.

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