FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSIS PROCEDURES.

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Intended use:

Lyophilised FPA controls, at a defined FibrinoPeptide A (FPA) concentration, as measured with the ZYMUTEST-FPA kit (ref RK016A), for the quality control of FPA antigenic assays. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

Reagents:

Each kit contains:

- 3 vials of 1 ml of FPA control I, lyophilised, at a high FPA concentration (CI)
- 3 vials of 1 ml of FPA control II, lyophilised, at a low FPA concentration (CII).

The FPA concentrations and the acceptance ranges for controls I and II, can present variations from lot to lot, but are indicated on the flyer provided with each set of controls.

Reagent Stability: In their original package, and before any use, when stored at 2-8 °C, the reagents are stable until the expiration date printed on the kit.

Excipients: Glycine, Sodium Chloride, stabilizers, Bovine Serum Albumine (BSA).

Viral safety: FPA controls and BSA were tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

Reconstitution and stability:

- **FPA Control I (high):** restore the vial with 1 ml sample diluent; let for 15 min at room temperature; mix gently until complete dissolution of the content.
- **FPA Control II (low):** restore the vial with 1 ml sample diluent; let for 15 min at room temperature; mix gently until complete dissolution of the content.

These controls are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota : FPA for FPA controls is obtained from human fibrinogen, prepared from plasmas which were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Preparation and results analysis using Zymutest FPA assay (RK016A):

No bentonite treatment is required for these controls. To **1 ml** of each reconstituted control, add exactly **0.1 ml of affinity purified** rabbit antibodies specific for FPA. Incubate for 1 hour at 37 °C. Controls are then tested exactly as samples using Zymutest FPA assay. From the curve obtained, deduce **directly** the FPA concentration in controls I and II.

Application

The FPA concentration, indicated on the flyer, for Control I and Control II, has been measured with the ZYMUTEST FPA (RK016A), Elisa kit.

The target values and the acceptance ranges, indicated on the flyer, are valid for ZYMUTEST FPA (ref RK016A) only. The results obtained can vary, when these controls are used with other kits according to the assay reactivity and its standardisation. These controls can nevertheless be used with different kits for quality control and testing run to run variability, but the acceptable values must be determined by each laboratory.

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