**DDIMER CONTROL PLASMA LOW**

FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSIS PROCEDURES.

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**Intended use**

Lyophilised plasma, at a defined DDimer concentration, as measured with commercially available methods and especially the VIDAS-DDimer assay

**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**Reagent**

Vials of 1 ml of Positive DDimer Control Human Plasma, lyophilised, at a low DDimer concentration (**DDimer control low**).

The DDimer concentration and the acceptance range, for this positive control, can present variations from lot to lot, but are exactly indicated on the flyer provided with each control.

The DDimer concentration measured with the Vidas DDimer kits is of about **0.6 µg/ml**.

**Reagent Stability:**

In its original package, and before any use, when stored at 2-8°C, the reagent is stable until the expiration date printed on the kit.

**Excipients:**

Glycine, stabilizers.

**Preparation:**

**DDimer Control low:** restore the vial with 1 ml of distilled water; mix with vortex during 5 seconds or until complete dissolution of the content, let for 15 min at room temperature;

This positive DDimer control plasma is stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in its original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for few min).

**Nota :** Plasmas used for preparing the DDimer control plasma were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Application**

The DDimer concentration, indicated on the flyer, for Positive DDimer Control plasma, has been measured with the Vidas DDimer, Elisa (Biomerieux).

The target value and the acceptance range, indicated on the flyer, are valid for the assay used. The results obtained can vary, when this control is used with other kits, according to the assay reactivity and its standardisation. This control can nevertheless be used with different DDimer assays, for the quality control, and for testing run to run variability, but the acceptance value and range must be determined by each laboratory.

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