

# NEW APPROACH FOR DETECTION OF HEPARIN DEPENDENT ANTIBODIES AND RISK ASSESSMENT FOR HEPARIN INDUCED THROMBOCYTOPENIA


Amiral J, Peyrafitte M, Catala M, Vissac AM  
HYPHEN BioMed Research, 95000 Neuville sur Oise (France)

## Introduction

- This assay uses the potential of immobilized and biologically active heparin to focus and catch antibody-protein (mainly PF4)-heparin complexes. It then mimics the conditions occurring *in vivo* when heparin dependent antibodies are generated and can induce Heparin Induced Thrombocytopenia (HIT).
- A new assay for measuring heparin dependent antibodies, involved in the development of HIT was developed.
- Various presentations are proposed for the measurement of total antibodies (IgGAM), or for specifically measuring IgG isotypes, or for the total isotyping of IgG, IgA and IgM isotypes.

## Assay principle

- Heparin, immobilized onto a solid reactive surface (plate or other), but «functionally available»:
  - Captures chemokines present into the patient plasma/serum (or supplied exogenously as a platelet lysate), and then forms the reactive auto-antigen, which binds heparin dependent antibodies.
  - Can also bind «heparin-protein-antibody» complexes present in blood circulation.
- «Functionally available» heparin uses one of the following coating procedures:
  - Protamine sulfate complexed with a large excess of heparin.
  - Streptavidin complexed with biotinylated heparin.
  - Heparin chemically coupled to a high molecular weight molecule (natural or synthetic) or polymer.



## Assay Protocol

- Plate coated with «functionally available» heparin.
- 1:100 (or more) diluted plasma or serum ± lysate.
- Second antibody (peroxidase labeled):
  - Anti-IgGAM (Screening) → HIT risk assessment
  - Anti-IgG (IgG isotype only) → Confirmation of HIT
  - Anti-IgG, Anti-IgA and anti-IgM (total isotyping) → Research studies
- TMB/H<sub>2</sub>O<sub>2</sub> substrate and reaction stopped with Sulfuric Acid.
- Measurement of absorbance at 450 nm (A450).

## References

- Amiral J et al : Platelet factor 4 complexed to heparin is the target for antibodies generated in heparin induced thrombocytopenia. *Thromb haemost* 1992; 68: 95-6.
- Warkentin TE et al. Heparin induced thrombocytopenia in patient treated with low molecular weight heparin or unfractionated heparin. *N Eng J Med* 1995; 332:1330-5.
- Amiral J et al. Antibodies to macromolecular platelet factor 4-heparin complexes in heparin induced thrombocytopenia: a study of 44 cases. *Thromb Haemost* 1995; 73: 21-28.
- Amiral J et al. Presence of autoantibodies to interleukin-8 or neutrophil-activating peptide-2 in patients with heparin associated thrombocytopenia. *Blood* 1996; 78:78-449 (abstract).
- Gruel Y. Thrombopénie induite par les héparines manifestations cliniques et physiopathologie. *Presse Med.* 1998; 27: S7-S12.
- Elalamy Y et al. Diagnostique et gestion des thrombopénies induites par l'héparine. *Rev Mal Respir* 1999, 16 : 961-974.
- Warkentin TE et al. Testing for heparin-induced thrombocytopenia antibodies. *Transfus Med Rev* 2006; 20: 259-72.
- Greinacher A. Heparin induced thrombocytopenia: frequency and pathogenesis. *Pathophysiol Haemost Thromb* 2006; 35:37-45.

## Conclusions

- New highly sensitive and specific assay for the diagnosis of heparin dependent antibodies involved in HIT, easy to perform and cost effective, offering automation possibilities.
- Good correlation with platelet aggregation tests and measurement of anti-H-PF4 antibodies.
- Potentially sensitive to the various antigenic targets for heparin dependent antibodies (studies in progress).
- Possible measurement of circulating complexes «heparin-protein-antibody» and assay mimicking the heparin dependent antibody binding mechanisms occurring *in vivo*.
- Very «flexible» assay principle for all laboratory immunological studies on heparin dependent antibodies, which can cause HIT.

## Results

Patients: Citrated plasmas from:

- 60 normal individuals
- 37 patients with a clinically diagnosed HIT (platelet course kinetics, positive platelet aggregation tests at low but not at high heparin concentration, recovery of platelet count following heparin withdrawal).

Table 1 :

A450 in normals and patients with HIT

Specificity	A450	SD
NI Plasmas (N=60)	<0.10	0.03
HIT Plasmas (N=37)	≥ 1.00	Range: 1.02 to >3.00

## Assay performances

