

Technical File

LIAPHEN™ vWF:Ag


Product Ref. N°	A120206-RUO
Kit presentation	Ready to use Latex automated method, R1. Reactive Buffer; R2: Latex
Status	RUO

Intended use:

Latex automated Immunoturbidimetric assay for measuring of von Willebrand Factor Antigen (vWF:Ag) in human citrated plasma.

Assay principle:

Immunoturbidimetric method, based on antigen-antibody reaction. The absorbance change is directly proportional to the amount of vWF:Ag in the sample.



Key Characteristics & Advantages:

vWF:Ag in practice:

▪ RUO Applications:

- For use in vWF research laboratories
- vWF functional or quantitative deficiency leads to von Willebrand disease (vWD) characterized by abnormal vWF adhesion activity.
- vWF deficiencies can develop in association with other pathologies, thus constituting an acquired von Willebrand disease.
- When vascular endothelium is affected, the vWF concentration can be increased in relation to inflammatory processes.
- Recent studies showed that very high vWF concentrations correlate with organ failure, liver disease severity, and might predict in-hospital survival in acute-on-chronic liver failure patients

▪ Market requirements & challenges:

- ✓ Rapid and performing vWF:Ag automated assay, highly stable, polyvalent for use on major automated coagulation analyzers. An extended dynamic range is of high value for the identification of vWD (low range) as well as for high to very high vWF concentrations

Easy and fast Automation:

- Simple, rapid (<5min), fully automatable on any coagulation analyzer.
- Method applications: Sysmex CS-series, STA-R® series, ACL TOP® family and BCS® XP.

OPTIMIZE YOUR ASSAY!

Optimized consumption with extended stability compared to other assays

- 15 days onboard (CS-series and STA-R®)
- Long shelf life at 2-8°C in original packaging

Extended measurement range (3-1600% vWF:Ag) useful to explore various contexts

No hook effect up to 1600%
Antibodies specific of vWF:Ag


Excellent Precision with very low CVs:

- Within-run < 4.6 % (example, STA-R and CS-5100)
- Between-run < 3.6% (example, STA-R and CS-5100)

Specific reliable assay, without interferences related to usual sample content (Hemoglobin, bilirubin, triglycerides, heparins and DOACs)

Standardized method for optimal routine testing

- Improved lot-to-lot variability; High quality raw materials used
- Calibrators/controls traceable to WHO International Standard




LIAPHEN™ vWF:Ag

Product Ref. N°	A120206-RUO
Kit presentation	Ready to use Latex automated method, R1. Reactive Buffer; R2: Latex
Status	RUO

Assay Characteristics*:


Multiplatform, easy to use, highly stable, low CVs , no interferences

Calibration curve:




Measurement range: 3-1600%
(automatic redilution > 150%)
Detection Limit: LLOD 0.7%
No hook effect: up to 1600%

Method Comparison (STA-R® and CS):




Good correlation with other device:
N=142
 $y = 1.06x - 3.9$
 $r = 0.998$



Good correlation with other device:
N=73
 $y = 1.03x + 0.49$
 $r = 0.998$

Analyzers Comparison:



Good correlation between analyzers:
N=55
 $y = 1.07x - 5.68$
 $r = 0.996$

Interferences and specificity:

High Specificity

- Polyclonal antibodies fully specific for vWF
- Expected normal concentrations: 62-169%
- Specificity verified on vWF deficient plasma: vWF:Ag %: < LOD

No interference of anticoagulant drugs

- Bias % was globally < 5% on controls (2IU/mL Heparins; 400ng/mL DOACs : Apixaban, Rivaroxaban, Dabigatran)

No interference of usual interferents and Plasma content

- Bias % was globally < 5% on controls (60 mg/dL Bilirubin, 1000mg/dL Hb or Intralipids)

Precision (CS-5100 / STA-R®):

Intra-series (N=40 determinations)

- | | |
|--|---|
| <ul style="list-style-type: none"> CS-5100 CV% 2.2% (vWF:Ag= 103%) CV% 4.6% (vWF:Ag= 40%) | <ul style="list-style-type: none"> STA-R® CV% 2.2% (vWF:Ag= 100%) CV% 3.6% (vWF:Ag= 38%) |
|--|---|

Inter-series (n=30 determinations, 5 days, 10 runs)

- | | |
|--|---|
| <ul style="list-style-type: none"> CS-5100 CV% 2.2% (vWF:Ag= 103%) CV% 2.6% (vWF:Ag= 39%) | <ul style="list-style-type: none"> STA-R® CV% 2.6% (vWF:Ag= 101%) CV% 3.6% (vWF:Ag= 38%) |
|--|---|

Kit presentation:



R1: 4 x 5 mL
R2: 4 x 6 mL

Example of stability and lot-to-lot variability:

- 3 lots on CS-5100:
- Normal [99-103%]
- Low [38-40%]

Lot-to-lot
variability
on same
controls

Shelf life in
original
packaging

- 24 months at 2-8°C

After
opening:
1 month 2-
8°C or 2
weeks RT

Onboard
15 days on
CS-5100

- Bias % was <5% on controls

- Fresh vs stored:
 Δ % (OD/min at
 $\approx$ 100%) was < 5%
(on CS)

* The values cited for the assay characteristics of the assay represent example of typical results and are not to be regarded as specifications for the test

Related products:

Needed	Product name	Ref N°:
Calibrator / Controls	BIOPHEN™ Plasma Calibrator	A222101
	BIOPHEN™ Normal Control Plasma	A223201
	BIOPHEN™ Abnormal Control Plasma	A223301
Buffer	Imidazole Buffer	AAR021A/K/L

For information	Product name	Ref N°:
ELISA	ZYMUTEST™ vWF	ARK030A-RUO
ELISA	ZYMUTEST™ vWF:CBA	ARK038A-RUO
Aggregometry	Ristocetin	AAG004K-RUO
Aggregometry	Lyophilized platelets	AAG006A or AAG006K

7768 Service Center Drive, West Chester, OH 45069

info@aniara.com • www.aniara.com

For Research Use Only. Not for Use in Diagnostic Procedures

ANIARA
DIAGNOSTICA
S.O. | HYPHEN BIOMED | MEDICAGO | XENOMETRIX