

## Data sheet SPP conjugate

<b>Cat. No.:</b>	<b>CT357</b>
<b>Format:</b>	5-plate
<b>Specificity:</b>	SPP provides the conjugate for the U-CyTech ELISA systems.
<b>Application:</b>	Reconstitute the lyophilized contents by injecting 500 µl of sterile distilled water into the vial. Mix the solution gently for approximately 15 seconds and allow it to stand for 1 minute at room temperature. Avoid vigorous shaking. Dilute the reconstituted SPP 100x in PBS with 0.5% BSA and 0.05% Tween-20. Apply 100 µl/well. Cover the plate and incubate 1 hour at 37°C. The product should be used in conjunction with TMB substrate (CT367).
<b>Contents:</b>	Lyophilized product; sufficient for 480 ELISA determinations. Buffer prior to lyophilization: 0.5 ml PBS + 125 mM trehalose + 1% BSA. Preservative free.
<b>Sterility:</b>	Membrane filtered (0.2 µm).
<b>Packaging:</b>	10 ml clear glass bottle.
<b>Storage:</b>	The lyophilized SPP can be best stored at ≤ -20°C for a defined length of time (expiry date indicated on vial). After reconstitution, the reagent is stable for at least 2 months at 4°C when kept sterile. However, it is strongly recommended to divide the reconstituted SPP preparations into small aliquots for single use. These aliquots should be stored at ≤ -20°C. Under these conditions the reagents are stable for minimal one year.
<b>General:</b>	A white pancake-like powder. After reconstitution: a clear solution. Do not use if the solution is or has become turbid.
<b>Safety:</b>	The product is not classified as hazardous according to Regulation (EC) no. 1272/2008 and to Directive 1999/45/EC and its amendments. Consult the material safety data sheet (MSDS) for further information regarding hazards and safe handling practices. The MSDS can be found on our website <a href="http://www.ucytech.com">www.ucytech.com</a> .

