**BIOPHEN CS-05(88)**  
**tPA and broad spectrum**  
**Chromogenic substrate**  
Vial of 25 mg/flacon de 25 mg

Ref. 229091

**STORE AT**  
**2-8 °C**

**For research use only**  
**Not for use in diagnosis**  
**procedures.**

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**CONSERVER A**  
**2-8 °C**

**Origin:**

Synthetic generic Chromogenic substrate, highly purified and stabilized.

**Specificity:**

Recommended substrate for tPA presenting a broad spectrum reactivity with many serine esterases, and particularly thrombin and plasmin.

**Presentation:**

Vial containing 25 mg of BIOPHEN CS-05(88), lyophilised in presence of Mannitol as a bulking agent.

**Reconstitution:**

According to the research protocol used, the BIOPHEN CS-05(88) chromogenic substrate can be restored with variable volumes of distilled water; for example 5 mL can be used for a substrate concentration of 5 mg/mL, or 20 mL for a substrate concentration of 1.25 mg/mL.

**Peptide sequence:** H-D-Ile-Pro-Arg-pNa-2HCl

**Molarity:** 49.5 µMoles

**Purity grade:** > 95%.

**Molecular weight:** 504.6.

**Free pNA content (measured on substrate at 2.5g/L)**  $OD_{405nm} \leq 0.30$

**Use:**

All research studies and protocols where a source of chromogenic substrate for tPA is required. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

**Respective reactivities:**

Plasmin	Thrombin	FXa	Urokinase	tPA
100	65	60	8	20

Assay conditions must be established for making the substrate totally specific for tPA.

**Stability:**

- Restored :**
- 7 days at room temperature
  - 3 months at 2-8 °C
  - **Do not freeze.**

**Lyophilized :** Until the expiration date printed on the vial.

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